

Supporting Information

Contents

1. Phantom Construction	3
2. MR Parameter Array Tables	4
3. Recommended Imaging Protocols	6
4. Materials Database	13
5. Fiducial Array/ Geometric Distortion Measurements	17
6. NMR Calibration Measurements	19
7. Field, Temperature, and Time Stability	23
8. Material Stability.....	27

Supporting Information Figures

FIGURE S1 A, The hemispherically-machined surface on the post matches with the inclined notches, B, to precisely locate the plates. Photos of blue fiducial, C, and red MnCl_2 array, D, spheres showing the sealing and mounting methods. D, inset shows micro-computed tomography image of contrast sphere.	3
FIGURE S2 Schematic of the fiducial array analysis showing input 3D image of the system phantom, cropped image with just the fiducial spheres, image of fiducial sphere, synthetic k-space image and real-space image used as a convolution mask, slice of 3D convolution image, and convolution profiles with fits used to obtain sphere center.	18
FIGURE S3 Fiducial analysis on a 1.5 T scanner with a gradient echo sequence. A, Coronal slice, B, axial slice, C, sagittal slice. D, magnified image of fiducial sphere image with ROI location after automated location. E, Normalized integrated intensity for all 57 fiducial spheres. F, Difference between center of mass and convolutional sphere centers.....	18
FIGURE S4 A, B, C, Geometric distortion of the 57 fiducial sphere centers in x, y and z directions, respectively.	19
FIGURE S5 NMR signal for Ni-12 versus inversion time for the T1-IR protocol along with fits to the model described in the text. Data for three consecutive measurements are shown along with the residuals for each measurement (top plot). The errors listed for T_1 and the inversion efficiency are the standard deviation of the 3 values obtained for each of the measurements.	21
FIGURE S6 NMR signal from Ni-12 using a CPMG sequence as a function of acquisition time along with exponential fits. Data for three consecutive measurements are shown along with the residuals for each measurement (top plot). The errors listed for T_2 are the standard deviation of the 3 values obtained for each of the measurements.	22
FIGURE S7 Magnetic field dependence of T_1 and T_2 for the MnCl_2 and NiCl_2 arrays.....	23
FIGURE S8 Temperature dependence of T_1 , T_2 for the Ni-12 solution measured in a metrology NMR at 3.0 T. The plot shows data from a flame-sealed borosilicate-capillary library sample over the course of 4 years.	24
FIGURE S9 Variation of normalized relaxation times with temperature for the NiCl_2 array at 3.0 T.	25
FIGURE S10 Temperature coefficient of spin relaxation times for the MnCl_2 array at 3.0 T.	26
FIGURE S11 Variation of relaxation times with temperature for the CuSO_4 fiducial solution at 3.0 T.....	26
FIGURE S12 Water mass uptake for various plastics: nylon/polyamide (PA), poly(methyl methacrylate) (PMMA), polycarbonate (PC), polyvinyl chloride (PVC), polyvinylidene fluoride (PVDF), polyphenylene sulfide (PPS), polypropylene (PP). The samples were 25 mm diameter, 6 mm thick disks.	27
FIGURE S13 Geometric distortion during water soaking of the same samples used in Figure S12. The horizontal line indicates the threshold for maintaining the specified geometric distortion of the phantom plates.	28

Supporting Information tables

TABLE S1: NiCl ₂ Array.....	4
TABLE S2: MnCl ₂ Array.....	4
TABLE S3: Proton Density Array.....	5
TABLE S4: Isotropic Volume Series	6
TABLE S5: Section Thickness Series	7
TABLE S6: Resolution Inset Series.....	8
TABLE S7: Proton Density and Signal to Noise Series	9
TABLE S8: T_1 Inversion Recovery Series.....	10
TABLE S9: T_1 Variable Flip Angle Series	11
TABLE S10: T_2 Series	12
TABLE S11: Starting Chemicals for Contrast Fluids.....	13
TABLE S12: NiCl ₂ Solutions.....	14
TABLE S13: MnCl ₂ Solutions.....	15
TABLE S14: Proton Density Solutions	16

1. Phantom Construction

Figure S1 shows photos of the plate locating mechanism required to achieve 0.1 mm accuracy of fiducial sphere centers. The hemispherically-machined surface on the post, A, matches with the inclined notches, B, to precisely locate the plates. Photos of blue fiducial, C, and red MnCl_2 , D, spheres showing the sealing and mounting methods. The fiducial sphere has a machined 10 mm internal diameter (ID) with a precise mounting and locating geometry. The base hemisphere and mounting stud are one machined piece, while the top hemisphere is glued into place, registered with a locating collar. The commercial polypropylene spheres used in the MR parameter arrays, do not have a precise spherical ID, as seen in the X-ray computed tomography image in the inset in D.




FIGURE S1 A, The hemispherically-machined surface on the post matches with the inclined notches, B, to precisely locate the plates. Photos of blue fiducial, C, and red MnCl_2 array, D, spheres showing the sealing and mounting methods. D, inset shows micro-computed tomography image of contrast sphere.

2. MR Parameter Array Tables

TABLE S1: NiCl₂ Array

Contrast ID	X (mm)	Y (mm)	Z (mm)	positional accuracy (mm)	ICP measured NiCl ₂ Concentration (mM)	NMR Measured T_1 @ 1.5 T (ms)	T_1 Standard Deviation (ms)	NMR Measured T_2 @ 1.5 T (ms)	T_2 Standard Deviation (ms)	NMR Measured T_1 @ 3.0 T (ms)	T_1 Standard Deviation (ms)	NMR Measured T_2 @ 3.0 T (ms)	T_2 Standard Deviation (ms)
Ni-1	0.0	56.5	50.0	±0.4	0.268	2 033	4.6	1 669	0.5	1 989	1.0	1 465	1.0
Ni -2	29.4	56.5	40.5	±0.4	0.556	1 489	1.4	1 244	0.6	1 454	2.5	1 077	1.8
Ni -3	47.6	56.5	15.5	±0.4	1.053	1 012	0.2	859.3	0.17	984.1	0.33	717.9	1.12
Ni -4	47.6	56.5	-15.5	±0.4	1.683	730.8	1.10	628.5	0.13	706	1.5	510.1	1.36
Ni -5	29.4	56.5	-40.5	±0.4	2.579	514.1	0.06	446.3	0.11	496.7	0.41	359.6	0.22
Ni -6	0.0	56.5	-50.0	±0.4	3.835	367.9	0.66	321.2	0.30	351.5	0.91	255.5	0.07
Ni -7	-29.4	56.5	-40.5	±0.4	5.674	260.1	0.04	227.7	0.07	247.13	0.086	180.8	0.04
Ni -8	-47.6	56.5	-15.5	±0.4	8.142	184.6	0.02	161.9	0.06	175.3	0.11	127.3	0.14
Ni -9	-47.6	56.5	15.5	±0.4	11.716	132.7	0.02	117.1	0.03	125.9	0.33	90.3	0.14
Ni -10	-29.4	56.5	40.5	±0.4	16.83	92.7	0.09	81.9	0.02	89.0	0.17	64.3	0.05
Ni-11	-20.0	56.5	20.0	±0.4	23.827	65.4	0.10	57.7	0.02	62.7	0.13	45.7	0.12
Ni-12	20.0	56.5	20.0	±0.4	33.922	46.32	0.010	41.0	0.01	44.53	0.090	31.86	0.02
Ni13	20.0	56.5	-20.0	±0.4	48.609	32.45	0.012	28.7	0.03	30.84	0.016	22.38	0.02
Ni-14	-20.0	56.5	-20.0	±0.4	68.884	22.859	0.043 7	20.2	0.01	21.719	0.0054	15.83	0.03

TABLE S2: MnCl₂ Array

Contrast ID	X (mm)	Y (mm)	Z (mm)	positional accuracy (mm)	ICP measured MnCl ₂ Concentration (mM)	NMR Measured T_1 @ 1.5 T (ms)	T_1 Standard Deviation (ms)	NMR Measured T_2 @ 1.5 T (ms)	T_2 Standard Deviation (ms)	NMR Measured T_1 @ 3.0 T (ms)	T_1 Standard Deviation (ms)	NMR Measured T_2 @ 3.0 T (ms)	T_2 Standard Deviation (ms)
Mn-1	0.0	16.5	50.0	±0.4	0.012	2 376	2.2	939.4	1.11	2 480	10.8	581.3	0.39
Mn -2	29.4	16.5	40.5	±0.4	0.019	2 183	1.2	594.3	0.32	2 173	14.7	403.5	0.55
Mn -3	47.6	16.5	15.5	±0.4	0.028	1 870	6.0	416.5	1.13	1 907	10.3	278.1	0.28
Mn -4	47.6	16.5	-15.5	±0.4	0.042	1 539	3.8	267.0	0.11	1 604	7.2	190.94	0.011
Mn -5	29.4	16.5	-40.5	±0.4	0.062	1 237	0.4	184.9	0.11	1 332	0.8	133.27	0.073
Mn -6	0.0	16.5	-50.0	±0.4	0.091	1 030	1.7	140.6	0.05	1 044	3.2	96.89	0.049
Mn -7	-29.4	16.5	-40.5	±0.4	0.131	752.2	1.12	91.76	0.029	801.7	1.70	64.07	0.034
Mn -8	-47.6	16.5	-15.5	±0.4	0.187	550.2	0.18	64.84	0.029	608.6	1.03	46.42	0.014
Mn -9	-47.6	16.5	15.5	±0.4	0.267	413.4	0.29	45.28	0.029	458.4	0.33	31.97	0.083
Mn -10	-29.4	16.5	40.5	±0.4	0.381	292.2	0.15	30.62	0.014	336.5	0.18	22.56	0.012
Mn -11	-20.0	16.5	20.0	±0.4	0.54	194.9	0.08	19.76	0.017	244.2	0.09	15.813	0.006 1
Mn -12	20.0	16.5	20.0	±0.4	0.762	160.2	0.23	15.99	0.012	176.6	0.09	11.237	0.005 7
Mn -13	20.0	16.5	-20.0	±0.4	1.08	106.4	0.02	10.47	0.006	126.9	0.03	7.911	0.003 7
Mn-14	-20.0	16.5	-20.0	±0.4	1.527	83.33	0.10	8.15	0.011	90.9	0.05	5.592	0.005 5

TABLE S3: Proton Density Array

Contrast ID	X (mm)	Y (mm)	Z (mm)	positional accuracy (mm)	D ₂ O %	NMR Measured T_1 @ 1.5 T (ms)	T_1 Standard Deviation (ms)	NMR Measured T_2 @ 1.5 T (ms)	T_2 Standard Deviation (ms)	NMR Measured T_1 @ 3.0 T (ms)	T_1 Standard Deviation (ms)	NMR Measured T_2 @ 3.0 T (ms)	T_2 Standard Deviation (ms)
PD-1	0.0	-23.5	50.0	±0.4	95	378.08	0.42	320.61	0.99	341.74	0.18	236.25	0.51
PD-2	29.4	-23.5	40.5	±0.4	90	371.76	0.10	311.21	0.52	336.69	0.08	233.59	0.27
PD-3	47.6	-23.5	15.5	±0.4	85	367.96	0.17	313.80	0.43	335.10	0.02	233.89	0.16
PD-4	47.6	-23.5	-15.5	±0.4	80	364.15	0.09	308.79	0.03	334.51	0.03	233.02	0.20
PD-5	29.4	-23.5	-40.5	±0.4	75	353.17	0.81	299.65	0.04	323.94	0.01	227.16	0.08
PD-6	0.0	-23.5	-50.0	±0.4	70	354.84	0.09	300.41	0.13	323.03	0.01	226.04	0.07
PD-7	-29.4	-23.5	-40.5	±0.4	65	365.43	0.18	310.05	0.28	343.58	0.05	240.48	0.12
PD-8	-47.6	-23.5	-15.5	±0.4	60	345.95	0.07	295.20	0.16	324.53	0.03	228.05	0.10
PD-9	-47.6	-23.5	15.5	±0.4	50	343.76	0.41	287.02	0.79	319.87	0.05	225.25	0.06
PD-10	-29.4	-23.5	40.5	±0.4	40	339.28	0.36	285.67	0.52	320.14	0.03	226.31	0.10
PD-11	-20.0	-23.5	20.0	±0.4	30	340.05	0.01	292.96	0.15	316.07	0.08	224.48	0.03
PD-12	20.0	-23.5	20.0	±0.4	20	328.78	0.08	284.15	0.06	309.31	0.02	220.68	0.09
PD-13	20.0	-23.5	-20.0	±0.4	10	329.23	0.05	287.47	0.05	311.68	0.03	224.09	0.06
PD-14	-20.0	-23.5	-20.0	±0.4	0	331.16	0.33	320.61	0.99	305.08	0.04	219.06	0.06

3. Recommended Imaging Protocols

TABLE S4: Isotropic Volume Series

Isotropic Volume Series	GE	Philips	Siemens
Sequence	3D/SPGR	3D/SPGR	3D/RF spoiled GRE
Scan Plane	COR	COR	COR
Scan Options	EDR	3D FFE; Fast = none	
Section Thickness / Gap (mm)	1	1	1
TR (ms)	5.9	10.0	6.3
TE (ms)	1.4 (minimum)	4.0	1.9
TI Values (ms)			
Flip Angle (deg)	10	10	10
ETL	1	1	1
Number of Averages	1	1	1
Matrix (FE)	256	256	256
Matrix (PE)	256	256	256
Matrix (SE) / # of Slices	256	256	256
Pixel Bandwidth (Hz)	488	434	650
Bandwidth (kHz) - GE	62.5		
FOV (FE, mm)	250	250	250
FOV (PE, mm)	250	250	250
Pixel Size (mm x mm)	0.98 x 0.98	0.98 x 0.98 x 0.98	0.98 x 0.98
PE Direction	RL	RL	RL
Notes	Resolutions are important for each of these protocols. Keep resolution and adjust TE/TR as necessary.	Resolutions are important for each of these protocols. Keep resolution and adjust TE/TR as necessary.	Resolutions are important for each of these protocols. Keep resolution and adjust TE/TR as necessary.
	Repeat for each plane (cor, sag, ax).	Repeat for each plane (cor, sag, ax).	Repeat for each plane (cor, sag, ax).
Series	Isotropic Volume	Isotropic Volume	Isotropic Volume
Approx. Acq. Time per Setting (min)	6.4	11.0	6.53
# Settings	1	1	1
Total Time for this Series (min)	6.4	11.0	6.53

TABLE S5: Section Thickness Series

Section Thickness Series	GE	Philips	Siemens
Sequence	2D/FSE	2D/SE	2D/TSE
Scan Plane	COR	COR	COR
Scan Options	EDR	MS SE; Fast=TSE	
Section Thickness / Gap (mm)	3/3, 5/5	3/3, 5/5	3/3, 5/5
TR (ms)	5000	2000	5000
TE (ms)	60	36	63
TI Values (ms)			
Flip Angle (deg)		90	
ETL	8	8	8
Number of Averages	2	1	2
Matrix (FE)	512	512	512
Matrix (PE)	256	256	256
Matrix (SE) / # of Slices	1	1	1
Pixel Bandwidth (Hz)	244	394	184
Bandwidth (kHz) - GE	62.5		
FOV (FE, mm)	250	250	250
FOV (PE, mm)	250	250	250
Pixel Size (mm x mm)	0.49 x 0.98	0.49 x 0.98	0.49 x 0.98
PE Direction	RL	RL	RL
Notes	Two series: one at 3mm slice and one at 5mm slice	Two series: one at 3mm slice and one at 5mm slice	Two series: one at 3mm slice and one at 5mm slice
Series	Section Thickness	Section Thickness	Section Thickness
Approx. Acq. Time per Setting (min)	5.3	1.5	2.47
# Settings	2	2	2
Total Time for this Series (min)	10.6	3.0	4.94

TABLE S6: Resolution Inset Series

Resolution Inset Series	GE	Philips	Siemens
Sequence	2D/FSE	2D/SE	2D/TSE
Scan Plane	COR	COR	COR
Scan Options	EDR	MS SE; Fast=TSE	
Section Thickness / Gap (mm)	4	4	4
TR (ms)	5000	5000	5000
TE (ms)	60	60	63
TI Values (ms)			
Flip Angle (deg)		90	
ETL	8	8	8
Number of Averages	2	2	2
Matrix (FE)	512	512	512
Matrix (PE)	512	504	512
Matrix (SE) / # of Slices	1	1	1
Pixel Bandwidth (Hz)	244	394	227
Bandwidth (kHz) - GE	62.5		
FOV (FE, mm)	250	250	250
FOV (PE, mm)	250	250	250
Pixel Size (mm x mm)	0.49 x 0.49	0.49 x 0.5	0.49 x 0.49
PE Direction	RL	RL	RL
Notes	Acquire 3 sections in one acquisition, spaced appropriately to cover the three resolution inserts (if extra coffins exist). Acquire both plate and coffin insets.	Acquire 3 sections in one acquisition, spaced appropriately to cover the three resolution inserts (if extra coffins exist). Acquire both plate and coffin insets.	Acquire 3 sections in one acquisition, spaced appropriately to cover the three resolution inserts (if extra coffins exist). Acquire both plate and coffin insets.
Series	Resolution Inset	Resolution Inset	Resolution Inset
Approx. Acq. Time per Setting (min)	10.7	10.6	5.27
# Settings	2	2	2
Total Time for this Series (min)	21.4	21.2	10.54

TABLE S7: Proton Density and Signal to Noise Series

PD & SNR Series	GE	Philips	Siemens
Sequence	2D/SE	2D/SE	2D/SE
Scan Plane	COR	COR	COR
Scan Options	EDR	MS SE; Fast=none	
Section Thickness / Gap (mm)	6	6	6
TR (ms)	5000	5000	5000
TE (ms)	10	10	10
TI Values (ms)			
Flip Angle (deg)		90	
ETL	1	1	1
Number of Averages	1	1	1
Matrix (FE)	256	256	256
Matrix (PE)	192	192	192
Matrix (SE) / # of Slices	1	1	1
Pixel Bandwidth (Hz)	244	213	153
Bandwidth (kHz) - GE	31.2		
FOV (FE, mm)	250	250	250
FOV (PE, mm)	200 (0.8 PFOV)	250	250
Pixel Size (mm x mm)	0.98 x 0.98	0.98 x 1.3	0.98 x 0.98
PE Direction	RL	RL	RL
Notes	Acquire twice and use NEMA MS 1-2008 Method 1 for SNR measurement.	Acquire as single-series, but two dynamics.	Acquire twice and use NEMA MS 1-2008 Method 1 for SNR measurement.
Series	PD & SNR	PD & SNR	PD & SNR
Approx. Acq. Time per Setting (min)	13.8	16.0	16.07
# Settings	2	2	2
Total Time for this Series (min)	27.7	32.0	32.14

TABLE S8: T_1 Inversion Recovery Series

T1 - VTI Series	GE	Philips	Siemens
Sequence	2D/FSE-IR	2D/IR-SK	2D/TSE-IR
Scan Plane	COR	COR	COR
Scan Options	EDR	2D IR; Fast=TSE	
Section Thickness / Gap (mm)	6	6	6
TR (ms)	4500	4500	4500
TE (ms)	Min Full (7.6)	7	6.9
TI Values (ms)	50, 75, 100, 125, 150, 250, 1000, 2000, 3000	35, 75, 100, 125, 150, 250, 1000, 1500, 2000, 3000	35, 75, 100, 125, 150, 250, 1000, 1500, 2000, 3000
Flip Angle (deg)			
ETL	3	6	6
Number of Averages	1	1	1
Matrix (FE)	256	256	256
Matrix (PE)	192	252	192
Matrix (SE) / # of Slices	1	1	1
Pixel Bandwidth (Hz)	391	436	279
Bandwidth (kHz) - GE	50		
FOV (FE, mm)	250	250	250
FOV (PE, mm)	200 (0.8 PFOV)	250	250
Pixel Size (mm x mm)	0.98 x 0.98	0.98 x 0.98	0.98 x 0.98
PE Direction	RL	RL	RL
Notes	Min TI allowed on GE: 50 ms Autoprescan with TI = 50 ms	Reconstruct Magnitude, Real, Imaginary images. Do "FullPrep" for each of these ten series.	Use SOS (sum-of-squares) multi-channel coil reconstruction, not ACC (adaptive coil combine). Known "line" artifacts result when using the ACC method.
Series	T1 VTI	T1 VTI	T1 VTI
Approx. Acq. Time per Setting (min)	4.0	3.5	3.02
# Settings	9	10	10
Total Time for this Series (min)	36.0	35.0	30.20

TABLE S9: T_1 Variable Flip Angle Series

T1 - VFA Series	GE	Philips	Siemens
Sequence	3D/FSPGR	3D/SPGR	3D/RF spoiled GRE
Scan Plane	COR	COR	COR
Scan Options	EDR/Z2	3D FFE; Fast=none	
Section Thickness / Gap (mm)	3/0	6/0	6/0 (or 3/0 if error)
TR (ms)	Min (6.0)	6.6	6.6
TE (ms)	Min (1.4)	1.8	2.44
TI Values (ms)			
Flip Angle (deg)	2, 5, 10, 20, 25, 30	2, 5, 10, 20, 25, 30	2, 5, 10, 20, 25, 30
ETL	1	1	1
Number of Averages	4	4	4
Matrix (FE)	256	256	256
Matrix (PE)	192	192	192
Matrix (SE) / # of Slices	34	28	32
Pixel Bandwidth (Hz)	488	904	280
Bandwidth (kHz) - GE	62.5		
FOV (FE, mm)	250	250	250
FOV (PE, mm)	250	250	250
Pixel Size (mm x mm)	0.98 x 0.98	0.98 x 1.3 x 6	0.98 x 0.98
PE Direction	RL	RL	RL
User CVs	Turbo = 0		
Notes	1) Yields 30 3 mm sections. 2) Autoprescan with FA=15. 3) Ensure gain settings do not vary between series, i.e., use Manual Prescan. 4) Include fiducial spheres above & below T1 spheres in scan volume. 5) Using Turbo=0 should maintain the same TE/TR for each FA.	1) Reconstructed at 3 mm. 2) Ensure gain settings do not vary between series, i.e., use MPS. Autoprescan using 15 deg FA. 3) Include fiducial spheres above & below T1 spheres in scan volume.	1) Reconstructed at 3 mm. 2) Ensure gain settings do not vary between series, i.e., use MPS. Autoprescan using 15 deg FA. 3) Include fiducial spheres above & below T1 spheres in scan volume.
Series	T1 VFA	T1 VFA	T1 VFA
Approx. Acq. Time per Setting (min)	2.6	3.0	1.23
# Settings	7	7	7
Total Time for this Series (min)	18.2	21.0	8.61

TABLE S10: T_2 Series

T2 Series	GE	Philips	Siemens
Sequence	2D/SE	2D/SE	2D/se_mc
Scan Plane	COR	COR	COR
Scan Options	EDR	2D SE; Fast=none	
Section Thickness / Gap (mm)	6	6	6
TR (ms)	5000	5000	5000
TE (ms)	15, 30, 45, 60; 25, 50, 75, 100; 40, 80, 120, 160	11 ms x 16 echoes	10 to 320 ms by 10 ms
TI Values (ms)			
Flip Angle (deg)			
ETL	1	1	1
Number of Averages	1	1	1
Matrix (FE)	256	256	256
Matrix (PE)	192	192	192
Matrix (SE) / # of Slices	1	1	1
Pixel Bandwidth (Hz)	156	172	227
Bandwidth (kHz) - GE	20.0		
FOV (FE, mm)	250	250	250
FOV (PE, mm)	200 (0.8 PFOV)	250	250
Pixel Size (mm x mm)	0.98 x 0.98	0.98 x 1.3	0.98 x 0.98
PE Direction	RL	RL	RL
Notes	Autoprescan with TE = 15 ms. Ensure gain settings do not vary between series, i.e., use Manual Prescan.	Ensure gain settings do not vary between series, i.e., use MPS.	
Series	T2	T1 VFA	T1 VFA
Approx Acq Time per Setting (min)	13.8	16.0	16.1
# Settings	3	1	1
Total Time for this Series (min)	41.5	16.0	16.1

4. Materials Database

TABLE S11: Starting Chemicals for Contrast Fluids

Chemical	Purity	Molecular Weight (g/mol)
CuSO ₄ •5H ₂ O	Technical grade	249.68
NiCl ₂ •6H ₂ O	99.999% trace metal basis	237.69
MnCl ₂ •4H ₂ O	99.99% trace metal basis	197.91
D ₂ O	D, 99.8%; chemical purity: 99.5%	20.03
H ₂ O	Deionized	18.02
Erioglaucine disodium salt	unknown	792.85
Tartrazine	≥85% Dye content	534.36
Allura Red AC	80% Dye content	496.42

Density of H₂O at 20 °C: 0.9982 g/cm³

Density of D₂O at 20 °C: 1.105 g/cm³

Fiducial solution [CuSO₄]:

Fiducial solutions were made at 20 °C by dissolving crystalline CuSO₄•5H₂O (Millipore-Sigma, St. Louis, MO, USA; Part number: 209198) into deionized water with color doped by erioglaucine disodium salt (FD&C Blue No. 1; Millipore-Sigma, St. Louis, MO, USA; Part number: 861146). The target concentration for Cu was 3.21 mM with a doping of 1.47 μM erioglaucine disodium salt. For making one liter of fiducial solution:

$$\text{Cu: } \frac{0.00321 \text{ mol}}{L} * 1 L * \frac{249.68 \text{ g}}{\text{mol}} = 0.8015 \text{ g}$$

$$\text{Erioglaucine disodium salt: } \frac{0.00000147 \text{ mol}}{L} * 1 L * \frac{792.85 \text{ g}}{\text{mol}} = 0.0012 \text{ g}$$

NiCl₂ Array:

NiCl₂ contrast array solutions were made at 20 °C from a stock solution of NiCl₂•6H₂O (Millipore-Sigma, St. Louis, MO, USA; Part number: 203866) with a concentration of 101 mM Ni⁺⁺, also made at 20 °C. Each contrast solution was color doped with erioglaucine disodium salt

and tartrazine (FD&C Yellow No. 5; Millipore-Sigma, St. Louis, MO, USA; Part number: T0388).

For making one liter of stock NiCl_2 solution:

$$\text{Ni: } \frac{0.101 \text{ mol}}{L} * 1 L * \frac{237.69 \text{ g}}{\text{mol}} = 24.0067 \text{ g}$$

To color dope each dilution of stock NiCl_2 solution, the following amount of dye was needed per dilution (assuming 100 mL volume):

$$\text{Erioglaurine disodium salt: } \frac{0.000035 \text{ mol}}{L} * 0.1 L * \frac{792.85 \text{ g}}{\text{mol}} = 0.00277 \text{ g}$$

$$\text{Tartrazine: } \frac{0.0000524 \text{ mol}}{L} * 0.1 L * \frac{534.36 \text{ g}}{\text{mol}} = 0.0028 \text{ g}$$

Based on the stock solution concentration, the following volumes are necessary for each of the dilutions (based on $C_{\text{initial}}V_{\text{initial}} = C_{\text{final}}V_{\text{final}}$), assuming 100 mL final volume:

TABLE S12: NiCl_2 Solutions

NiCl ₂ Contrast Number	Desired Concentration (mM)	Volume of stock solution needed (mL)
Ni-1	0.299	0.27
Ni 1-2	0.62	0.56
Ni -3	1.074	0.98
Ni -4	1.716	1.56
Ni -5	2.624	2.39
Ni -6	3.908	3.56
Ni -7	5.724	5.21
Ni -8	8.293	7.55
Ni -9	11.925	10.85
Ni -10	17.062	15.53
Ni -11	24.327	22.14
Ni -12	34.6	31.48
Ni -13	49.129	44.70
Ni -14	69.677	64.40

MnCl₂ Contrast Array:

MnCl₂ contrast array solutions were made at 20 °C from a stock solution of MnCl₂•4H₂O (Millipore-Sigma, St. Louis, MO, USA; Part number: 221279) with a concentration of 25 mM, also made at 20 °C. Each contrast solution was color doped with Allura Red AC (FD&C Red 40; Millipore-Sigma, St. Louis, MO, USA; Part number: 458848). For making one liter of stock MnCl₂ solution:

$$\text{Mn: } \frac{0.025 \text{ mol}}{L} * 1 L * \frac{197.91 \text{ g}}{\text{mol}} = 4.9478 \text{ g}$$

To color dope each dilution of stock MnCl₂ solution, the following amount of dye was needed per dilution (assuming 100 mL volume):

$$\text{Allura Red AC: } \frac{0.0001 \text{ mol}}{L} * 0.1 L * \frac{496.42 \text{ g}}{\text{mol}} = 0.00496 \text{ g}$$

Based on the stock solution concentration, the following volumes are necessary for each of the dilutions (based on $C_{\text{initial}}V_{\text{initial}} = C_{\text{final}}V_{\text{final}}$), assuming 100 mL final volume:

TABLE S13: MnCl₂ Solutions

MnCl ₂ Contrast Number	Desired Concentration (mM)	Volume of stock solution needed (mL)
Mn-1	0.01206	0.0482
Mn -2	0.01989	0.0796
Mn -3	0.03076	0.1230
Mn -4	0.04663	0.1865
Mn -5	0.06879	0.2752
Mn -6	0.1001	0.4004
Mn -7	0.14442	0.5777
Mn -8	0.20705	0.8282
Mn -9	0.29565	1.1826
Mn -10	0.42094	1.6838
Mn -11	0.59814	2.3926
Mn -12	0.84873	3.3949
Mn -13	1.20313	4.8125
Mn -14	1.70431	6.8172

Proton Density Contrast Array [D₂O/H₂O]:

Proton density (PD) contrast array solutions were made at 20 °C from solutions of D₂O, H₂O and stock solution of NiCl₂ (noted previously in *T*₁ contrast array section). Each proton density contrast solution was made independently, and color doped with tartrazine.

To color dope each proton density contrast solution, the following amount of dye was needed per dilution (assuming 100 mL volume):

$$\text{Tartrazine: } \frac{0.0001 \text{ mol}}{L} * 0.1 L * \frac{534.36 \text{ g}}{\text{mol}} = 0.0053 \text{ g}$$

To reduce *T*₁ of water for acceptable scan times, each proton density contrast solution is doped with 0.1 mL of NiCl₂ stock solution (101 mM) for a final concentration of (assuming 100 mL volume):

$$\text{Ni: } 0.1 \text{ mL} * 101 \text{ mM} = 100 \text{ mL} * x \text{ mM} \rightarrow x = 0.101 \text{ mM}$$

Because of the addition of water as part of the NiCl₂ stock solution, it is important to subtract that volume from the water used for the remainder of the proton density contrast solution.

Based on the stock solution concentration, the following volumes are necessary for each of the dilutions (based on $C_{\text{initial}}V_{\text{initial}} = C_{\text{final}}V_{\text{final}}$), assuming 100 mL final volume:

TABLE S14: Proton Density Solutions

PD Contrast Number	Desired H ₂ O/D ₂ O (%/%)	Volume of water after NiCl ₂ stock solution (mL)
PD-1	5/95	4.9
PD-2	10/90	9.9
PD-3	15/85	14.9
PD-4	20/80	19.9
PD-5	25/75	24.9
PD-6	30/70	29.9
PD-7	35/65	34.9
PD-8	40/60	39.9
PD-9	50/50	49.9
PD-10	60/40	59.9

PD-11	70/30	69.9
PD-12	80/20	79.9
PD-13	90/10	89.9
PD-14	100/0	99.9

5. Fiducial Array/ Geometric Distortion Measurements

Figure S2 schematically shows the fiducial array analysis process. The input 3D gradient echo image is cropped after a manual coarse alignment of the ROIs to the fiducial spheres. The 3D cropped image is then convolved with a synthetic fiducial sphere mask. A fiducial sphere image is shown as a comparison to the mask. Both real and synthetic images have partial volume effects and finite k-space sampling structure (Gibbs ringing). A slice of the 3D convolution is shown along with lines scans that are fit with a Gaussian model to determine sphere centers.

Figure S3 shows the 10 mm ROIs after sphere location and after the translation, rotation, and scaling, required to minimize $\sum_i (\vec{R}_{ai} - \vec{R}_{pi})^2$, is applied. Here, \vec{R}_{ai} is the apparent position and \vec{R}_{pi} is the prescribed center position of the i^{th} sphere. Once the sphere centers have been located, the integrated signal $I_i = \sum_j I_{ij}$, shown in Figure S3E, can be calculated as a measure of image uniformity. The sum of voxel intensities I_{ij} for the i^{th} sphere is over all voxels, j , whose centers are within 6 mm from the sphere center. For the scan shown in Figure S3, the anterior and posterior spheres have distinctly lower signal. The sphere centers can be calculated with a variety of algorithms. Figure S3F shows a comparison between sphere centers determined by fitting the convolution profiles and by calculating the center of mass $\vec{R}_{ai} = \sum_j \vec{R}_{aj} I_{ij} / I_i$, where the sum is over the same the voxel set used to calculate the integrated intensity. As seen in Figure 3F, the agreement is within 0.04 mm.

Figure S4 shows the distortion analysis from a scan that has non-uniform gradient corrections. In this case, they appear to be only done in the scan plane. The corrected x and z distortions (in-plane distortion) have standard deviations of 0.11 mm and 0.22 mm, respectively, whereas the y distortion perpendicular to the scan plane has a standard deviation of 0.63 mm. A completely independent analysis on the same data set, which located the sphere centers by finding their boundary region, gave standard deviations in the in-plane x, z directions of 0.17 mm, 0.17 mm, respectively, and the out-of-plane y direction of 0.66 mm.




FIGURE S2 Schematic of the fiducial array analysis showing input 3D image of the system phantom, cropped image with just the fiducial spheres, image of fiducial sphere, synthetic k-space image and real-space image used as a convolution mask, slice of 3D convolution image, and convolution profiles with fits used to obtain sphere center.




FIGURE S3 Fiducial analysis on a 1.5 T scanner with a gradient echo sequence. A, Coronal slice, B, axial slice, C, sagittal slice. D, magnified image of fiducial sphere image with ROI location after automated location. E, Normalized integrated intensity for all 57 fiducial spheres. F, Difference between center of mass and convolutional sphere centers.




FIGURE S4 A, B, C, Geometric distortion of the 57 fiducial sphere centers in x , y and z directions, respectively.

6. NMR Calibration Measurements

Here we give an overview of SI-traceable NMR measurement protocols of proton relaxation times, fully described in NIST SP-250-97 (1). Briefly, NMR samples were prepared by injecting the contrast fluids into Teflon or flame sealed borosilicate glass capillaries for single measurements and long-term monitoring, respectively. Complete libraries containing 3 flame sealed capillaries of all solutions were made. Temperature measurements were made using a fiberoptic thermometer

placed 15 mm from the center of the RF coils. The fiberoptic thermometer was calibrated against two NIST-calibrated platinum resistance thermometers. The uncertainty in the temperature is typically ± 0.3 C.

T_1 Measurement

The longitudinal relaxation time, T_1 , is defined as the exponential recovery time of the local magnetization along the applied field direction after it is driven out of equilibrium. NMR inversion recovery (IR) experiments were conducted by monitoring longitudinal relaxation of ^1H nuclei at a frequency corresponding to 128 MHz (± 0.1 MHz) and 64 MHz (± 0.1 MHz) at 20 °C. The NMR IR experiment utilized a composite 180° RF pulse for inversion to minimize errors relating to the calculation of the 90° and 180° times. Each reported T_1 value is obtained from a collection of 20 inversion time (TI) spectra, where TI is the time between the end of the 180° RF inversion pulse and the start of 90° RF pulse, which tips the spins back into the transverse plane for data acquisition. The TI values are chosen to approximately span the range from $0.01 T_1$ to $5 T_1$. Each inversion time spectrum is the Fourier transform of the average of eight phase-cycled free-induction decay measurements. The delay between measurements is required to be greater than five times the T_1 value, if it is determined that the delay between measurements is shorter than five times the T_1 value, the experiment is repeated with an appropriate length delay. The complex spectra were phase shifted so that the real component gave a symmetric Lorentzian peak. Peak integration was determined by fitting the spectrum with a Lorentzian and numerically integrating ± 10 linewidths. Integrated values from each inversion time spectra, shown in Figure S5, were fit with an exponential function of the form:

$$S(TI) = A \left(1 - (1 + \delta) e^{-\frac{TI}{T_1}} \right).$$

Where $S(TI)$ is the integrated signal intensity for an inversion time TI and δ is a fit parameter that indicates the degree of inversion in the peak integration data with an ideal value of $\delta = 1$. Each experiment was repeated three times without adjusting the sample, nor repeating the shim. Fitted T_1 times were then averaged together, and deviation is reported for the three measurements from the mean T_1 value (SD_{3R}). NMR details of the experiments can be found in (2) with composite pulse details found in (3).




FIGURE S5 NMR signal for Ni-12 versus inversion time for the T1-IR protocol along with fits to the model described in the text. Data for three consecutive measurements are shown along with the residuals for each measurement (top plot). The errors listed for T_1 and the inversion efficiency are the standard deviation of the 3 values obtained for each of the measurements.

T_2 Measurement

NMR Carr-Purcell-Meiboom-Gill (CPMG) experiments were conducted by monitoring transverse relaxation of ^1H nuclei at a frequency corresponding to 128 MHz (± 0.1 MHz) and 64 MHz (± 0.1 MHz) at 20 °C. The NMR experiment consists of a 90° RF pulse to tip the spins, about the x-axis, into the transverse plane and then a train of 180° RF refocusing pulses about the y-axis until the desired acquisition time, $t_a = n(2\tau_{cp} + t_{180})$, is reached. Here, n is the number of 180° RF refocusing pulses, $\tau_{cp} = 1$ ms is the delay before and after the 180° pulse, and t_{180} is the duration of the 180° pulse. Each reported T_2 value is obtained from a collection of 20 CPMG sequences where the number of refocusing pulses is increased to vary t_a . t_a values are logarithmically-spaced and are chose to have a maximum time that insures observation of at least three decades of signal decrease. For these measurements, the signal-to-noise ratio varied from 2000 to 5000. Each CPMG spectrum is the Fourier transform of the average of eight phase-cycled free-induction decay measurements. The delay between measurements is required to be greater than five times

the T_1 value (as determined by NMR-IR measurements), if it is determined that the delay between measurements is shorter than five times the T_1 value, the experiment is repeated with an appropriate length delay. The complex spectra were phase shifted so that the real component gave a symmetric Lorentzian peak. Peak integration was determined by fitting the spectrum with a Lorentzian and numerically integrating ± 10 linewidths. Integral values from each CPMG spectra $S(t_a)$ vs. acquisition time are shown in Figure S6 along with fits to an exponential function of the form:

$$S(t_a) = S_0 e^{-t_a/T_2}.$$

Each experiment was repeated three times without adjusting the sample, nor repeating the shim. Fitted T_2 times were then averaged together, and deviation is reported for the three measurements from the mean T_2 value (SD_{3R}). The transverse relaxation time, T_2 , will be dependent on the refocusing time, $2\tau_{cp} + t_{180}$, and, therefore, T_2 must be reported with the refocusing time used. NMR details of NMR-CPMG experiments can be found in (2) with original references (4-6).




FIGURE S6 NMR signal from Ni-12 using a CPMG sequence as a function of acquisition time along with exponential fits. Data for three consecutive measurements are shown along with the residuals for each measurement (top plot). The errors listed for T_2 are the standard deviation of the 3 values obtained for each of the measurements.

7. Field, Temperature, and Time Stability

The proton relaxation times, as measured by NMR, of the MR-parameter solutions exhibit considerable magnetic field and temperature dependence. The field dependencies of T_1 , T_2 , taken from the measured relaxation times given in Tables 1, 2, are shown in Figure S7. Here, the




FIGURE S7 Magnetic field dependence of T_1 and T_2 for the MnCl_2 and NiCl_2 arrays.

normalized change in proton relaxation time per unit field, in units of $\%/T$, is given by

$$\frac{100}{1.5T} \frac{T_{1,2}(3T) - T_{1,2}(1.5T)}{T_{1,2}(1.5T)}$$

and is only a coarse measurement of the field dependence. Note that the NiCl_2

T_1 times decrease with increasing field, while the MnCl_2 T_1 , along with most tissues, increase with field. Given the large field dependence, especially of T_2 , small variations in imaging field of commercial scanners operating at nominally the same field values, must be taken into account.

An example of the temperature dependence at 3.0 T and 4-year stability of the Ni-12 solution is shown in Figure S8. The variations in T_1 , T_2 over the typical MRI bore temperature range (18 C to 26 C) are 1.3 % and 1.6 %, respectively. The repeatability in the measurement of the same sample for 3 measurements made 49 months apart is less than 0.1 %.




FIGURE S8 Temperature dependence of T_1 , T_2 for the Ni-12 solution measured in a metrology NMR at 3.0 T. The plot shows data from a flame-sealed borosilicate-capillary library sample over the course of 4 years.

The temperature variation of the proton spin relaxation rates in the NiCl_2 array are complex, as seen in Figure S9, exhibiting a minimum around 20 C for the higher concentrations, while exhibiting a more linear behavior for lower concentrations. The minimum in T_1 is characteristic of the resonance frequency, in this case 128 MHz, being aligned with a characteristic fluctuation rate. The minimum in relaxation times is useful to provide less temperature variation, but it is more difficult to provide thermal correction formulas since the form of the temperature variation changes with concentration and field.




FIGURE S9 Variation of normalized relaxation times with temperature for the NiCl_2 array at 3.0 T.

The MnCl_2 solutions exhibit a linear variation of T_1 and T_2 over the measurement range of 16 C to 26 C, with both T_1 and T_2 increasing with temperature. T_1 and T_2 have temperature coefficients of $\sim 3\%/C$ and $2\%/C$, respectively, at 3.0 T as shown in Figure S10. The CuSO_4 fiducial solution exhibits a similar linear dependence of T_1 and T_2 on temperature, as seen in Figure S11, with temperature coefficients of $2.72\%/C$ and $2.65\%/C$, respectively




FIGURE S10 Temperature coefficient of spin relaxation times for the MnCl₂ array at 3.0 T.




FIGURE S11 Variation of relaxation times with temperature for the CuSO₄ fiducial solution at 3.0 T.

8. Material Stability

Properties of materials used in the phantom were studied to determine water uptake and dimensional stability. A plot of water uptake is shown in Figure S12 and dimensional changes during water absorption are shown in Figure S13 for a variety of plastics. Water uptake can lead to undesirable water loss/exchange with the surrounding environment and to geometric distortion. Based on this data, PVC and PP were selected for the fiducial and contrast arrays, respectively, and PPS was selected for the components requiring high dimensional tolerance (the structural plates). Note, since there are no standard properties for many plastics, this data is for guidance only. Each specific material used in the phantom must be tested.




FIGURE S12 Water mass uptake for various plastics: nylon/polyamide (PA), poly(methyl methacrylate) (PMMA), polycarbonate (PC), polyvinyl chloride (PVC), polyvinylidene fluoride (PVDF), polyphenylene sulfide (PPS), polypropylene (PP). The samples were 25 mm diameter, 6 mm thick disks.




FIGURE S13 Geometric distortion during water soaking of the same samples used in Figure S12. The horizontal line indicates the threshold for maintaining the specified geometric distortion of the phantom plates.

REFERENCES

1. Boss MA, Dienstfrey AM, Gimbutas Z, Keenan KE, Splett JD, Stupic KF, Russek SE. Magnetic Resonance Imaging Biomarker Calibration Service: Proton Spin Relaxation Times. Volume Special Publication (NIST SP) - 250-97: NIST; 2018.
2. Berger S, Braun S. 200 and More NMR Experiments: A Practical Course: Wiley-VCH; 2004. 854 p.
3. Levitt MH, Freeman R. Nmr Population-Inversion Using A Composite Pulse. J Magn Reson 1979;33(2):473-476.
4. Hahn EL. SPIN ECHOES. Physical Review 1950;80(4):580-594.
5. Carr HY, Purcell EM. Effects Of Diffusion On Free Precession In Nuclear Magnetic Resonance Experiments. Physical Review 1954;94(3):630-638.
6. Meiboom S, Gill D. Modified Spin-Echo Method For Measuring Nuclear Relaxation Times. Review of Scientific Instruments 1958;29(8):688-691.